

MTT Assay

Day 1: Change media. For confluent cell lines, reduce cell number by half.

Day 2: Plate cells in 96 well plates

1. trypsinize and pellet cells
2. resuspend pellet in 3 or 5 ml of 10% serum media depending on pellet size
3. syringe 5-10 times through 18 gauge needle to get single cell suspension, vortex
4. add 10 μ l cell suspension to 90 μ l trypan blue, count with hemacytometer
 $\# \text{ of cells/ml} = \# \text{ of cells}/4 \times \text{dilution factor} \times 10^4$
dilution factor=10
5. dilute cells to $1 \times 10^5/\text{ml}$ in 10% serum media, vortex
6. count cells again (no trypan blue- undiluted)
7. dilute cells to desired concentration, vortex
for 4,000 cells/well, dilute cells to $4 \times 10^4/\text{ml}$
8. plate 100 μ l per well (10 rows of 6, leave empty wells around the outside)
always vortex before counting or diluting to assure even cell suspension

Day 3: Add drug

1. plan dilution scheme based on desired concentration of drug
2. dilute drug in serum-free media, **vortex in between each serial dilution**
3. add 100 μ l per well
4. add 200 μ l sterile water per well to each empty well

Day 8: Add MTT

1. make solution of 2 mg MTT/ml serum free media (heated to 37° C)- need 3.3 ml/plate
2. sterilize MTT solution with syringe filter
3. add 50 μ l per well, including a blank row
4. leave in incubator for 4 hours
5. suck off liquid using aspirator in the lab- make sure not to suck off purple solid in the bottom of the wells
6. add 100 μ l solubilizing solution, pipet up and down to mix
for 400 μ l solubilizing solution: 92 μ l DI water, add 8 μ l HCl to the water, add that to 300 μ l propanol
7. read plate at 490 with reference wavelength of 650, automix once before reading

10 mL = ~~2.3~~ 2.3 mL H₂O
100 μ l HCl 0.1N
7.5 mL iso

20 mL 4.6
200 10N
15 mL